

PTIP antibody

Catalog Number: 114306

Product name

PTIP antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PTIP Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 100-105 kDa

Preparation

This antibody was obtained by immunization of PTIP recombinant protein (Accession Number: BC033781). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114306(PAXIP1 antibody) at dilution of 1:500



IP Result of anti-PAXIP1 (IP:Catalog No:114306, 4ug; Detection:Catalog No:114306 1:500) with HeLa cells lysate 2400ug.