

FOXO4 antibody

Catalog Number: 110706

Product name

FOXO4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FOXO4 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, mouse uterus tissue. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 65kd

Preparation

This antibody was obtained by immunization of FOXO4 recombinant protein (Accession Number: NM_005938). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110706(FOXO4 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human kidney using Catalog No:110706(FOXO4 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:110706(FOXO4 antibody) at dilution of 1:100 (under 40x lens)