

# Pentraxin 3 antibody



Catalog Number: 113706

## Product name

Pentraxin 3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Pentraxin 3 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in HepG2 cells, Jurkat cells. Observed molecular weight by Western-blot: 42 kDa

## Preparation

This antibody was obtained by immunization of Pentraxin 3 recombinant protein (Accession Number: NM\_002852). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113706 (PTX3 antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells using Catalog No:113706 (PTX3 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)