

Anti-Cdk4 Rabbit antibody



Catalog Number: 166706

Product name

Anti-Cdk4 Rabbit antibody

Specificity

Human

Antibody description

Cdk4 Rabbit polyclonal antibody

Preparation

Antigen: Peptide sequence around aa.295~299(L-H-K-D-E) derived from Human CDK4.

Formulation

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, IHC-P


Dilutions

WB:1/1000IHC:1/50-1/100

Validations



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue



using CDK4 Antibody .

Western blot analysis of extracts from HeLa, 293T and HUVEC cells using CDK4 rabbit pAb (1:1000 diluted). Predicted band size: 34KDa. Observed band size: 34KDa.