

INPP5D Polyclonal Antibody



Catalog Number: 160016

Product name

INPP5D Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to INPP5D

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 959-1188 of human INPP5D (NP_005532.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations



Western blot - INPP5D Polyclonal Antibody

Western blot analysis of extracts of Jurkat cells, using INPP5D antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - INPP5D Polyclonal Antibody

Immunofluorescence analysis of NIH/3T3 cells using INPP5D antibody at dilution of 1:100. Blue: DAPI for nuclear staining.