

General Information

Gene Name Synonym

Collagenase 2; Matrix metalloproteinase-8

Protein Construction

A DNA sequence encoding the pro form of mouse MMP8 (NP_032637.3) (Phe 21-Ser 465) was expressed and purified, with an initial Met at the N-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant mouse MMP8 consisting of 446 amino acids and has a calculated molecular mass of 52 kDa. It migrates as an approximately 65 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse MMP-8 / CLG1 Protein SDS-PAGE