

Product name

NAGK antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

NAGK Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse testis tissue. Observed molecular weight by Western-blot: 37kd

Preparation

This antibody was obtained by immunization of NAGK recombinant protein (Accession Number: NM_017567). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113016(NAGK antibody) at dilution of 1:500