

PLVAP antibody



Catalog Number: 114016

Product name

PLVAP antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PLVAP Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, human kidney tissue, L02 cells, mouse liver tissue. Observed molecular weight by Western-blot: 50-60kd

Preparation

This antibody was obtained by immunization of PLVAP recombinant protein (Accession Number: NM_031310). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114016(PLVAP antibody) at dilution of 1:800