

SULF2 Polyclonal Antibody



Catalog Number: 166016

Product name

SULF2 Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to SULF2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 726-870 of human SULF2 (NP_061325.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - SULF2 Polyclonal Antibody

Western blot analysis of extracts of mouse brain, using SULF2 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.