

ZNF649 antibody



Catalog Number: 117016

Product name

ZNF649 antibody

Immunogen

[Human ZNF649 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

ZNF649 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, mouse heart tissue. Observed molecular weight by Western-blot: 50-57 kDa

Preparation

This antibody was obtained by immunization of ZNF649 recombinant protein (Accession Number: NM_023074). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:117016(ZNF649 Antibody) at dilution of 1:300