

Product name

Amyloid beta (A4) precursor protein antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Amyloid beta (A4) precursor protein Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 102 kDa

Preparation

This antibody was obtained by immunization of Amyloid beta (A4) precursor protein recombinant protein (Accession Number: BC004369). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108016(APP antibody) at dilution of 1:400