

Anti-ERK1/2 Rabbit antibody



Catalog Number: 169016

Product name

Anti-ERK1/2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

ERK1/2 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human ERK2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


1:2000-1:10000

1:50

1:100

1:20

Validations



Western blot detection of p42 MAPK (ERK2) in K562, Rat Brain, C6, 3T3, HeLa cell lysates using p42 MAPK (ERK2) Rabbit mAb(1:1000 diluted). Predicted band size: 43,41 kDa. Observed band size: 44,42kDa.