

# CD27 antibody



Catalog Number: 109016

## Product name

CD27 antibody

## Specificity

Human; other species not tested.

## Antibody description

CD27 Rabbit Polyclonal antibody. Positive FC detected in Raji cells. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 50 kDa

## Preparation

This antibody was obtained by immunization of CD27 recombinant protein (Accession Number: NM\_001242). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, FC

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:109016(CD27 antibody) at dilution of 1:300



1X10<sup>6</sup> Raji cells were stained with 0.2ug CD27 antibody (Catalog No:109016, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.