

# ORAI1 antibody



Catalog Number: 113416

## Product name

ORAI1 antibody

## Specificity

Human, Mouse, Rat, Zebrafish; other species not tested.

## Antibody description

ORAI1 Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells. Positive WB detected in Jurkat cells, A375 cells, MCF7 cells, mouse ovary tissue, mouse spleen tissue, PC-3 cells. Observed molecular weight by Western-blot: 45-47 kDa

## Preparation

This antibody was obtained by immunization of ORAI1 recombinant protein (Accession Number: NM\_032790). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113416(ORAI1 antibody) at dilution of 1:800



Immunofluorescent analysis of HEK-293 cells using Catalog No:113416(ORAI1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)