

# Anti-NRG2 antibody



Catalog Number: 176416

## Product name

WB, IHC-P

Anti-NRG2 antibody

## Dilutions

## Specificity

WB: 1:500-1:2000

Human, Mouse, Rat

IHC-P: 1:100-300

## Antibody description

Rabbit polyclonal antibody to NRG2

## Validations

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human nrg2. aa range:361-410

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal


## Ig Type

Rabbit IgG

## Applications



Fig1:: Western Blot analysis of HT29 cells using Neuregulin-2 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at



1:20000

Fig2:: Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:100