

# Anti-Glutathione S Transferase mu 1 (1H4) Mouse antibody



Catalog Number: 168416

## Product name

Mouse IgG1

Anti-Glutathione S Transferase mu 1 (1H4) Mouse antibody

## Applications

WB, IHC-P, FC

## Specificity

Human, Monkey

## Dilutions

WB: 1/500 - 1/2000

## Antibody description

Glutathione S Transferase mu 1 (1H4) Mouse monoclonal antibody

IHC: 1/200 - 1/1000

FC: 1/200 - 1/400

## Preparation

Antigen: Purified recombinant fragment of human GSTM1 expressed in E. Coli.

ELISA: 1/10000

## Validations

## Formulation

Ascitic fluid containing 0.03% sodium azide.


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of Glutathione S Transferase mu in MCF7, Jurkat and HeLa cell lysates using Glutathione S Transferase mu mouse mAb (dilution 1:1000). Predicted band size: 26 Kda. Observed band size: 26KDa.