

IP6K2 antibody



Catalog Number: 111816

Product name

IP6K2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

IP6K2 Rabbit Polyclonal antibody. Positive WB detected in human adrenal gland tissue, human heart tissue, K-562 cells. Observed molecular weight by Western-blot: 49kd

Preparation

This antibody was obtained by immunization of IP6K2 recombinant protein (Accession Number: XM_017006585). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human adrenal gland tissue were subjected to SDS PAGE followed by western blot with Catalog No:111816(IP6K2 antibody) at dilution of 1:500