

ABHD6 antibody



Catalog Number: 107816

Product name

ABHD6 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ABHD6 Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells. Observed molecular weight by Western-blot: 35-40 kDa

Preparation

This antibody was obtained by immunization of ABHD6 recombinant protein (Accession Number: NM_001320126). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107816(ABHD6 Antibody) at dilution of 1:300