

Human PEPD/PRD recombinant protein



Catalog Number: 519816

General Information

Gene Name Synonym

Xaa-Pro dipeptidase, Imidodipeptidase, PeptidaseD, Prolinedipeptidase, PRD, PEPD.

Protein Construction

A DNA sequence encoding the Human PEPD/PRD Ala2-Lys493 is expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM PB \square 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human PEPD/PRD recombinant protein