

DUSP13 antibody

Catalog Number: 110116

Product name

DUSP13 antibody

Immunogen

[Human DUSP13 Recombinant protein \(GST tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

DUSP13 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer. Positive WB detected in human skeletal muscle tissue, mouse testis tissue. Observed molecular weight by Western-blot: 22kd,36kd

Preparation

This antibody was obtained by immunization of DUSP13 recombinant protein (Accession Number: NM_001320843). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:110116(DUSP13 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:110116(DUSP13 antibody) at dilution of 1:100 (under 25x lens)