

ARPC1A Polyclonal Antibody



Catalog Number: 162116

Product name

ARPC1A Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to ARPC1A

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 171-350 of human ARPC1A (NP_006400.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100

Validations



Western blot - ARPC1A Polyclonal Antibody

Western blot analysis of extracts of U-251MG cells, using ARPC1A antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.