

# KPNA6 antibody

Catalog Number: 112116

## Product name

KPNA6 antibody

## Immunogen

[Human KPNA6 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

KPNA6 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells, HeLa cells. Positive IHC detected in human prostate cancer tissue. Observed molecular weight by Western-blot: 55-60 kDa

## Preparation

This antibody was obtained by immunization of KPNA6 recombinant protein (Accession Number: NM\_012316). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:112116(KPNA6 antibody) at dilution of 1:50



(under 10x lens)

BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:112116(KPNA6 Antibody) at dilution of 1:600