

# PPP4C antibody

Catalog Number: 114116

## Product name

PPP4C antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PPP4C Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, Jurkat cells, mouse kidney tissue. Positive IP detected in mouse kidney tissue. Observed molecular weight by Western-blot: 35 kDa

## Preparation

This antibody was obtained by immunization of PPP4C recombinant protein (Accession Number: NM\_002720). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114116(PPP4C antibody) at dilution of 1:1000



IP Result of anti-PPP4C (IP:Catalog No:114116, 4ug; Detection:Catalog No:114116 1:1000) with mouse kidney tissue lysate 4000ug.