

# RBBP7 Polyclonal Antibody



Catalog Number: 165116

## Product name

RBBP7 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to RBBP7

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human RBBP7 (NP\_002884.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IF, IP, ChIP

## Dilutions


WB 1:500 - 1:2000

IF 1:50 - 1:100

IP 1:50 - 1:200


ChIP 1:50 - 1:200

## Validations



Western blot - RBBP7 Polyclonal Antibody


Western blot analysis of extracts of various cell lines, using RBBP7 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 1s.

Immunofluorescence - RBBP7 Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using RBBP7 antibody . Blue: DAPI for nuclear staining.




Immunoprecipitation - RBBP7 Polyclonal Antibody

# RBBP7 Polyclonal Antibody



Catalog Number: 165116

Immunoprecipitation analysis of 200ug extracts of MCF7 cells using 1ug RBBP7 antibody . Western blot was performed from the immunoprecipitate using RBBP7 antibody at a dilution of 1:1000.



Chromatin Immunoprecipitation - RBBP7 Polyclonal Antibody

Chromatin immunoprecipitation analysis of extracts of 293T cells, using RBBP7 antibody and rabbit IgG. The amount of immunoprecipitated DNA was checked by quantitative PCR. Histogram was constructed by the ratios of the immunoprecipitated DNA to the input.