

# Anti-fezf2 antibody

Catalog Number: 175116

---

## Product name

Anti-fezf2 antibody

## Specificity

Zebrafish

## Antibody description

Rabbit polyclonal antibody to fezf2

## Preparation

This antigen of this antibody was recombinant protein within human fezf2 aa 211-385.

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 50% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB: 1:500-1:1000

## Validations




Fig1: Western blot analysis of FEZF2 on Zebrafish tissue lysates using anti-FEZF2 antibody at 1/500 dilution.