

THYN1 antibody



Catalog Number: 116116

Product name

THYN1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

THYN1 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in HEK-293 cells, mouse lung tissue. Observed molecular weight by Western-blot: 28-32kd

Preparation

This antibody was obtained by immunization of THYN1 recombinant protein (Accession Number: NM_199297). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116116(THYN1 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using Catalog No:116116 and Rhodamine-labeled goat anti-rabbit IgG (red).