

Anti-DFNA5/GSDME Rabbit antibody



Catalog Number: 167116

Product name

Anti-DFNA5/GSDME Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB, IP

Antibody description

DFNA5/GSDME Rabbit polyclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Recombinant protein of human DFNA5/GSDME

IP: 1/20

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Validations



Western blot detection of DFNA5/GSDME in A549 cell lysates using DFNA5/GSDME Rabbit pAb(1:1000 diluted). Predicted band size:55KDa.Observed band size:55KDa.