

AP1G2 antibody



Catalog Number: 108116

Product name

AP1G2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AP1G2 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, MCF-7 cells, mouse brain tissue. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of AP1G2 recombinant protein (Accession Number: XM_005268180). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:108116(AP1G2 Antibody) at dilution of 1:300