

Anti-Prolactin/PRL Mouse antibody



Catalog Number: 169116

Product name

Anti-Prolactin/PRL Mouse antibody

Specificity

Human

Antibody description

Prolactin/PRL Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human PRL (AA: 65-173) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG2b


Applications

ICC/IF, FC


Validations



Western blot analysis using PRL mAb against HEK293 (1) and PRL (AA)





Immunofluorescence analysis of MCF-7 cells using



PRL mouse mAb (green). Blue

Western blot analysis using PRL mAb against



Flow cytometric analysis of HeLa cells using PRL mouse mAb (green) and negative control (red).

Immunofluorescence analysis of HeLa cells using PRL mouse mAb (green). Blue