

Anti-Prolactin/PRL Mouse antibody



Catalog Number: 169116

Product name

Anti-Prolactin/PRL Mouse antibody

Specificity

Human

Antibody description

Prolactin/PRL Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human PRL (AA: 65-173) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG2b


Applications

ICC/IF, FC


Validations



Western blot analysis using PRL mAb against HEK293 (1) and PRL (AA



Immunofluorescence analysis of MCF-7 cells using




PRL mouse mAb (green). Blue

Western blot analysis using PRL mAb against


Anti-Prolactin/PRL Mouse antibody



Catalog Number: 169116



Flow cytometric analysis of Hela cells using PRL mouse mAb (green) and negative control (red).



Immunofluorescence analysis of Hela cells using PRL mouse mAb (green). Blue