

CDCA5 antibody

Catalog Number: 109116

Product name

CDCA5 antibody

Specificity

Human; other species not tested.

Antibody description

CDCA5 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of CDCA5 recombinant protein (Accession Number: NM_080668). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109116(CDCA5 Antibody) at dilution of 1:1000