

# JMJD8 Polyclonal Antibody



Catalog Number: 160516

## Product name

Rabbit IgG

JMJD8 Polyclonal Antibody

## Applications

## Specificity

WB

Human, Mouse, Rat

## Dilutions

## Antibody description

WB 1:500 - 1:2000

Polyclonal antibody to JMJD8

## Validations

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-285 of human JMJD8 (NP\_001005920.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type



## Western blot - JMJD8 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using JMJD8 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3s.