

Human CD123/IL3RA (His & Fc Tag) recombinant protein



Catalog Number: 500516

General Information

Protein Construction

A DNA sequence encoding the human IL3RA (NP_002174.1) extracellular domain (Met 1-Arg 305) was fused with the C-terminal polyhistidine-tagged Fc region of human IgG1 at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Activity

Measured by its binding ability in a functional ELISA. Immobilized human IL3 at 20 µg/ml (100 µl/well) can bind human IL3RA with a linear ranger of 0.064-1.6 µg/ml.

Purity

> 75 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Thr 19

Molecular Mass

The recombinant human IL3RA/Fc is a disulfide-linked homodimer. The reduced monomer comprises 535 amino acids and has a predicted molecular mass of 61 kDa. As a result of glycosylation, the apparent molecular mass of rh IL3RA/Fc monomer is approximately 80-90 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human IL3RA / CD123 Protein (His & Fc Tag)