

SLC24A5 Polyclonal Antibody



Catalog Number: 162516

Product name

SLC24A5 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to SLC24A5

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 220-300 of human SLC24A5 (NP_995322.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - SLC24A5 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using SLC24A5 antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.