

# RAD51AP1 antibody



Catalog Number: 114516

## Product name

RAD51AP1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RAD51AP1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse liver tissue, rat liver tissue. Observed molecular weight by Western-blot: 38 kDa

## Preparation

This antibody was obtained by immunization of RAD51AP1 recombinant protein (Accession Number: NM\_001130862). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114516(RAD51AP1 antibody) at dilution of 1:300