

# UBE2B antibody



Catalog Number: 116516

## Product name

UBE2B antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

UBE2B Rabbit Polyclonal antibody. Positive WB detected in L02 cells, Transfected HEK-293 cells. Positive IHC detected in human ovary tumor tissue. Observed molecular weight by Western-blot: 17 kDa

## Preparation

This antibody was obtained by immunization of UBE2B recombinant protein (Accession Number: NM\_003337). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:10-1:100

## Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:116516(UBE2B Antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:116516(UBE2B antibody) at dilution of 1:10 (under 10x lens)