

Anti-A1BG Rabbit antibody



Catalog Number: 167516

Product name

Anti-A1BG Rabbit antibody

Specificity

Human

Antibody description

A1BG Rabbit polyclonal antibody

Preparation

Antigen: Synthesized peptide derived from A1BG .
at AA range: 100-180

Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA
and 0.02% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ELISA


Dilutions

WB: 1/500 - 1/2000


IHC-p: 1/100-300

ELISA: 1/40000


Validations



Western Blot analysis of extracts from K562 cells, using A1BG Polyclonal Antibody. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-uterus, antibody was diluted at




1:100

Immunohistochemical analysis of paraffin-embedded human-ovary, antibody was diluted at


Anti-A1BG Rabbit antibody

Catalog Number: 167516



1:100

Immunohistochemical analysis of paraffin-embedded human-uterus, antibody was diluted at



1:100

Immunohistochemical analysis of paraffin-embedded human-ovary, antibody was diluted at 1:100