

# S1 tag antibody



Catalog Number: 107516

---

## Product name

S1 tag antibody

## Specificity

Human; other species not tested.

## Antibody description

S1 tag Mouse Monoclonal antibody. Positive WB detected in Recombinant protein.

## Preparation

This antibody was obtained by immunization of Peptide. Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:5000-1:50000

## Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:117337(S1 tag Antibody) at dilution of 1:10000