

# MYCBPAP antibody



Catalog Number: 112916

## Product name

MYCBPAP antibody

## Specificity

Human; other species not tested.

## Antibody description

MYCBPAP Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 108 kDa

## Preparation

This antibody was obtained by immunization of MYCBPAP recombinant protein (Accession Number: NM\_032133). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:112916(MYCBPAP antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells, using MYCBPAP antibody Catalog No:112916 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).