

# Anti-PREB Rabbit antibody



Catalog Number: 168916

## Product name

Rabbit IgG

Anti-PREB Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF, IP

Human

## Dilutions

## Antibody description

1:1000-1:5000

PREB Rabbit monoclonal antibody

1:20-1:50

## Preparation

1:20-1:50


Antigen: Recombinant protein of human PREB

1:20

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Validations



## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of PREB in K562,C6,3T3,HeLa cell lysates using PREB Rabbit mAb(1:1000 diluted).Predicted band size:46kDa.Observed band size:46kDa.

## Clonality

Monoclonal

## Ig Type