

# DHDH antibody

Catalog Number: 109916

## Product name

DHDH antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

DHDH Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 39kd

## Preparation

This antibody was obtained by immunization of DHDH recombinant protein (Accession Number: NM\_014475). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109916(DHDH antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109916(DHDH antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109916(DHDH antibody) at dilution of 1:100 (under 40x lens)