

SESN3 antibody

Catalog Number: 115216

Product name

SESN3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SESN3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A431 cells, A549 cells, HEK-293 cells, HepG2 cells, mouse liver tissue, mouse ovary tissue. Observed molecular weight by Western-blot: 51-59 kDa

Preparation

This antibody was obtained by immunization of SESN3 recombinant protein (Accession Number: BC017296). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115216(SESN3 antibody) at dilution of 1:500