

# Anti-Transferrin (4H5) Mouse antibody



Catalog Number: 168216

## Product name

Anti-Transferrin (4H5) Mouse antibody

## Specificity

Human

## Antibody description

Transferrin (4H5) Mouse monoclonal antibody

## Preparation

Antigen: Transferrin human.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300.

## Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Mouse IgG1


## Applications

ELISA

## Validations



Standard Curve for Transferrin: Capture Antibody Mouse mAb (2D5-E9) to Transferrin at 4u03bcg/ml and Detector Antibody Mouse mAb (4H5-H7) was used for detecting.



Observed Transferrin levels in several Human serum samples at different dilution.