

CAV3 Monoclonal Antibody



Catalog Number: 160616

Product name

CAV3 Monoclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Monoclonal antibody to CAV3

Preparation

Antigen: Recombinant protein of human CAV3

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB 1:500 - 1:2000

IP 1:20 - 1:50

Validations



Western blot - CAV3 Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using CAV3 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.