

# PBX4 antibody



Catalog Number: 113616

---

## Product name

PBX4 antibody

## Immunogen

[Human PBX4 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

PBX4 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, A2780 cells, COLO 320 cells, Jurkat cells, SKOV-3 cells. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of PBX4 recombinant protein (Accession Number: NM\_025245). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113616(PBX4 Antibody) at dilution of 1:600