

UNC5A-Specific antibody



Catalog Number: 116616

Product name

UNC5A-Specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UNC5A-Specific Rabbit Polyclonal antibody. Positive FC detected in MCF-7 cells. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 103 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_133369). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116616(UNC5A-Specific antibody) at dilution



of 1:100

1X10⁶ MCF-7 cells were stained with 0.2ug UNC5A-Specific antibody (Catalog No:116616, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.