

# Anti-PPCS Rabbit antibody



Catalog Number: 167616

## Product name

Anti-PPCS Rabbit antibody

## Specificity

Human

## Antibody description

PPCS Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human PPCS

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions

WB: 1/1000

IP: 1/20

## Validations



Western blot detection of PPCS in HeLa, Jurkat cell lysates using PPCS Rabbit mAb(1:1000 diluted). Predicted band size: 34kDa. Observed band size: 34kDa.