

# Anti-TNS4/CTEN Rabbit antibody



Catalog Number: 168616

## Product name

Anti-TNS4/CTEN Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat

## Applications

WB

## Antibody description

TNS4/CTEN Rabbit monoclonal antibody

## Dilutions

1:1000

## Preparation

Antigen: A synthetic peptide of human TNS4

## Validations

## Formulation


Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal



Western blot detection of TNS4/CTEN in K562 cell lysates using TNS4/CTEN Rabbit mAb(1:1000 diluted). Predicted band size:77kDa. Observed band size:90kDa.