

C18orf55 antibody



Catalog Number: 108616

Product name

C18orf55 antibody

Specificity

Human; other species not tested.

Antibody description

C18orf55 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue, human heart tissue, human kidney tissue, human ovary tissue, human spleen tissue.

Preparation

This antibody was obtained by immunization of C18orf55 recombinant protein (Accession Number: NM_014177). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human liver using Catalog No:108616(C18orf55 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:108616(C18orf55 antibody) at dilution of 1:100 (under 40x lens)