

# GULP1 antibody

Catalog Number: 111316

## Product name

GULP1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

GULP1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas. Positive WB detected in HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 45kd

## Preparation

This antibody was obtained by immunization of GULP1 recombinant protein (Accession Number: BC001103). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111316(GULP1 antibody) at dilution of 1:1500



Immunohistochemical of paraffin-embedded human pancreas using Catalog No:111316(GULP1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human pancreas using Catalog No:111316(GULP1 antibody) at dilution of 1:50 (under 40x lens)